

Rat Disintegrin And Metalloproteinase Domain-Containing Protein 33 ELISA Kit | ADAM33 elisa kit

Rat Disintegrin And Metalloproteinase Domain-Containing Protein 33 ELISA Kit

CATALOG # AAA210399	REACTIVITY Rat
-------------------------------	--------------------------


PRODUCT SYNONYMS

Disintegrin And Metalloproteinase Domain-Containing Protein 33; Rat Disintegrin And Metalloproteinase Domain-Containing Protein 33 ELISA Kit; ADAM33 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rat
Sequence Length	812
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	187.5 pg/ml-3000 pg/ml
Sensitivity	10 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for ADAM33 elisa kit	Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of Disintegrin and metalloproteinase domain-containing protein 33 in the sample, This Disintegrin and metalloproteinase domain-containing protein 33 Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus Disintegrin and metalloproteinase domain-containing protein 33 concentration. The concentration of Disintegrin and metalloproteinase domain-containing protein 33 in the samples is then determined by comparing the O.D. of the samples to the standard curve.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	540344525
NCBI Accession #	NP_001269376.1
NCBI GenBank Nucleotide #	NP_001269376.1
UniProt Accession #	Q9BZ11
NCBI Official Full Name	disintegrin and metalloproteinase domain-containing protein 33 isoform gamma
UniProt Protein Name	Disintegrin and metalloproteinase domain-containing protein 33

UniProt Gene Name	ADAM33
UniProt Synonym Gene Names	C20orf153
UniProt Entry Name	ADA33_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/rat-elisa-kits/disintegrin-and-metalloproteinase-domain-containing-protein-33/210399>