

Rat Cytokine-Induced Neutrophil Chemoattractant-3 ELISA Kit

Rat Cytokine-Induced Neutrophil Chemoattractant-3 ELISA Kit

CATALOG #
AAA210627

REACTIVITY
Rat


PRODUCT SYNONYMS

Cytokine-Induced Neutrophil Chemoattractant-3; Rat Cytokine-Induced Neutrophil Chemoattractant-3 ELISA Kit; Cytokine-Induced Neutrophil Chemoattractant-3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Rat
Samples	Serum, plasma, cell culture supernate and other biological fluids
Assay Type	Sandwich
Detection Range	62.5 pg/ml-1000 pg/ml
Sensitivity	1 pg/ml

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Cytokine-Induced Neutrophil Chemoattractant-3 elisa kit

Intended Use: The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of CYTOKINE-INDUCED NEUTROPHIL CHEMOATTRACTANT-3 in the sample, This CYTOKINE-INDUCED NEUTROPHIL CHEMOATTRACTANT-3 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus CYTOKINE-INDUCED NEUTROPHIL CHEMOATTRACTANT-3 concentration. The concentration of CYTOKINE-INDUCED NEUTROPHIL CHEMOATTRACTANT-3 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.