

Mouse anti-Human CD117/c-Kit Monoclonal Antibody | anti-KIT antibody

CD117/c-Kit (Marker for Gastrointestinal Stromal Tumors)

CATALOG #	HOST	REACTIVITY	APPLICATIONS
AAA215052	Mouse	Human	ELISA

PRODUCT SYNONYMS

CD117/c-Kit; Monoclonal Antibody; CD117/c-Kit (Marker for Gastrointestinal Stromal Tumors); p145; Steel Factor Receptor; Stem Cell Factor Receptor (SCF-Receptor); Mast Cell Growth Factor Receptor; anti-KIT antibody

PRODUCT SPECIFICATIONS

Host	Mouse
Reactivity	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone Number	KIT/2669
Specificity	This MAb recognizes a protein of 145kDa, identified as CD117/p145kit. It is found on a wide variety of tumor cells including follicular and papillary carcinoma of thyroid, adenocarcinomas from endometrium, lung, ovary, pancreas, and breast as well as malignant melanoma, endodermal sinus tumor, and small cell carcinoma. However, anti-CD117 has been particularly useful in differentiating gastrointestinal stromal tumors from Kaposi's sarcoma, tumors of smooth muscle origin, fibromatosis, and neural tumors of the GI tract. Anti-CD117 is also useful in recognizing myeloblasts in bone marrow biopsy and clot section.
Purity/Purification	Purified Ab with BSA and Azide at 200ug/ml OR Purified Ab WITHOUT BSA and Azide at 1.0mg/ml
Form/Format	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.
Sequence Length	976
Applicable Applications for anti-KIT antibody	ELISA
Positive Control	Gastrointestinal Stromal Tumor (GIST) or testicular germ cell tumor. Melanocytes in the basal layer of the epidermis and mast cells in the dermis of normal skin.
Immunogen	Recombinant full-length human KIT protein
Cellular Localization	Cell Surface and Cytoplasmic
Hu-chromosome Location	4q12
Preparation and Storage	Antibody with azide: store at 2 to 8 degree C. Antibody without azide: store at -20 to -80 degree C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

IMAGES



Application Data

Analysis of Protein Array containing more than 19,000 full-length human proteins using CD117 Mouse Monoclonal Antibody (KIT/2669). Z- and S- Score: The Z-score represents the strength of a signal that a monoclonal antibody (MAb) (in combination with a fluorescently-tagged anti-IgG secondary antibody) produces when binding to a particular protein on the HuProt™ array. Z-scores are described in units of standard deviations (SD's) above the mean value of all signals generated on that array. If targets on HuProt™ are arranged in descending order of the Z-score, the S-score is the difference (also in units of SD's) between the Z-score. S-score therefore represents the relative target specificity of a MAb to its intended target. A MAb is considered to specific to its intended target, if the MAb has an S-score of at least 2.5. For example, if a MAb binds to protein X with a Z-score of 43 and to protein Y with a Z-score of 14, then the S-score for the binding of that MAb to protein X is equal to 29.



SDS-PAGE

SDS-PAGE Analysis Purified CD117 Mouse Monoclonal Antibody (KIT/2669). Confirmation of Purity and Integrity of Antibody.

ADDITIONAL INFORMATION

References	1. Harvey JM, et. al. Journal of Clinical Oncology 1999; 17(5), 1474-1481. 2. Zafrani B, et. al. Histopathology 2000; 37(6), 536-545.
------------	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	4557695
NCBI GeneID	3815
NCBI Accession #	NP_000213.1
NCBI GenBank Nucleotide #	NP_000213.1
UniProt Accession #	P10721
Molecular Weight	~145kDa
NCBI Official Full Name	mast/stem cell growth factor receptor Kit isoform 1
NCBI Official Synonym Full Names	KIT proto-oncogene, receptor tyrosine kinase
NCBI Official Symbol	KIT
NCBI Official Synonym Symbols	PBT; SCFR; C-Kit; CD117; MASTC
NCBI Protein Information	mast/stem cell growth factor receptor Kit
UniProt Protein Name	Mast/stem cell growth factor receptor Kit
UniProt Gene Name	KIT
UniProt Synonym Gene Names	SCFR; SCFR; PBT
UniProt Entry Name	KIT_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.