

Human NOS3/eNOS ELISA Kit | NOS3/eNOS elisa kit

Human NOS3/eNOS (Nitric Oxide Synthase 3, Endothelial) ELISA Kit

CATALOG #
AAA21676

REACTIVITY
Human


PRODUCT SYNONYMS

NOS3/eNOS; Human NOS3/eNOS (Nitric Oxide Synthase 3; Endothelial) ELISA Kit; NOS3/eNOS elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This kit recognizes Human NOS3/eNOS in samples. No significant cross-reactivity or interference between Human NOS3/eNOS and analogues was observed.
Sequence Length	1203
Assay Type	Quantitative Sandwich
Samples	Serum, plasma and other biological fluids
Detection Range	62.5-4000pg/mL
Sensitivity	37.5pg/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, mid range and high level Human NOS3/eNOS were tested 20 times on one plate, respectively.
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, mid range and high level Human NOS3/eNOS were tested on 3 different plates, 20 replicates in each plate.
Preparation and Storage	Store at 4 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for NOS3/eNOS elisa kit	Principle of the Assay: This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Human NOS3/eNOS. Standards or samples are added to the micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Human NOS3/eNOS and AvidinHorseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Human NOS3/eNOS, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm ± 2 nm. The OD value is proportional to the concentration of Human NOS3/eNOS. You can calculate the concentration of Human NOS3/eNOS in the samples by comparing the OD of the samples to the standard curve. Intended Uses: This ELISA kit applies to the in vitro quantitative determination of Human NOS3/eNOS concentrations in serum, plasma and other biological fluids.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	40254422
NCBI GeneID	4846
NCBI Accession #	NP_000594.2
NCBI GenBank Nucleotide #	NM_000603.4
UniProt Accession #	P29474; Q13662; Q14251; Q14434; Q548C1; Q6GSL5; Q9UDC6; A0S0A7; A0S0A8; A8KA63; B2RCQ1; E9PFR2
Molecular Weight	67,876 Da
NCBI Official Full Name	nitric oxide synthase, endothelial isoform 1
NCBI Official Synonym Full Names	nitric oxide synthase 3 (endothelial cell)
NCBI Official Symbol	NOS3
NCBI Official Synonym Symbols	eNOS; ECNOS
NCBI Protein Information	nitric oxide synthase, endothelial; EC-NOS; NOS type III; NOSIII; cNOS; constitutive NOS; endothelial NOS; nitric oxide synthase 3 transcript variant eNOS-delta20; nitric oxide synthase 3 transcript variant eNOS-delta20-21; nitric oxide synthase 3 transcript variant eNOS-delta21
UniProt Protein Name	Nitric oxide synthase, endothelial
UniProt Gene Name	NOS3
UniProt Synonym Gene Names	cNOS; eNOS; NOSIII
UniProt Entry Name	NOS3_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/nos3-enos/21676>