

Human Mitogen Activated Protein Kinase 12 (MAPK12) RTU ELISA Kit | MAPK12 rtu elisa kit

Human Mitogen Activated Protein Kinase 12 (MAPK12) ELISA Kit

CATALOG #
AAA216924

REACTIVITY
Human


PRODUCT SYNONYMS

Mitogen Activated Protein Kinase 12 (MAPK12); Human Mitogen Activated Protein Kinase 12 (MAPK12) ELISA Kit; ERK6; PRKM12; p38 Gamma; SAPK3; Stress-activated protein kinase 3; Extracellular signal-regulated kinase 6; Mitogen-activated protein kinase p38 gamma; MAPK12 rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	357
Assay Type	Sandwich
Samples	Tissue Homogenates or Other Biological Fluids.
Detection Range	0.156-10ng/mL
Sensitivity	0.072ng/mL
RTU ELISA Kit's Advantage	The Ready-To-Use (RTU) Kits require fewer steps in setup and take less time to run. Detection Solutions are pre-diluted and total incubation time is shorter than many other ELISA Kits.
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for MAPK12 rtu elisa kit

Test Principle: The microtiter plate provided in this kit has been pre-coated with an antibody specific to MAPK12. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to MAPK12. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain MAPK12, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of MAPK12 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	735997389
NCBI GeneID	6300
NCBI Accession #	NP_001290181.1
NCBI GenBank Nucleotide #	NM_001303252.2
UniProt Accession #	P53778; Q14260; Q6IC53; Q99588; Q99672
Molecular Weight	40,808 Da
NCBI Official Full Name	mitogen-activated protein kinase 12 isoform 2
NCBI Official Synonym Full Names	mitogen-activated protein kinase 12
NCBI Official Symbol	MAPK12
NCBI Official Synonym Symbols	ERK3; ERK6; ERK-6; SAPK3; PRKM12; SAPK-3; MAPK 12; P38GAMMA
NCBI Protein Information	mitogen-activated protein kinase 12
UniProt Protein Name	Mitogen-activated protein kinase 12
UniProt Gene Name	MAPK12
UniProt Synonym Gene Names	ERK6; SAPK3; MAP kinase 12; MAPK 12; ERK-6; MAP kinase p38 gamma

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-rtu-elisa-kits/mitogen-activated-protein-kinase-12-mapk12/216924>