

Human Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4 (NOX4) RTU ELISA Kit | NOX4 rtu elisa kit

Human Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4 (NOX4) ELISA Kit

CATALOG #
AAA217426

REACTIVITY
Human


PRODUCT SYNONYMS

Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4 (NOX4); Human Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4 (NOX4) ELISA Kit; KOX; KOX-1; RENOX; NADPH Oxidase 4; NOX4 rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	538
Assay Type	Sandwich
Samples	Tissue Homogenates, Cell Lysates or Other Biological Fluids.
Detection Range	0.156-10ng/mL
Sensitivity	0.051ng/mL
RTU ELISA Kit's Advantage	The Ready-To-Use (RTU) Kits require fewer steps in setup and take less time to run. Detection Solutions are pre-diluted and total incubation time is shorter than many other ELISA Kits.
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for NOX4 rtu elisa kit

Test Principle: The microtiter plate provided in this kit has been pre-coated with an antibody specific to NOX4. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to NOX4. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain NOX4, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of NOX4 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	219842346
NCBI GeneID	50507
NCBI Accession #	NP_001137308.1
NCBI GenBank Nucleotide #	NM_001143836.2
UniProt Accession #	Q9NPH5; Q5K3R4; Q5K3R5; Q5K3R6; Q5K3R8; Q7Z7G3; Q86V92; A8K715; B7Z520; E7EMD7
Molecular Weight	59,915 Da
NCBI Official Full Name	NADPH oxidase 4 isoform b
NCBI Official Synonym Full Names	NADPH oxidase 4
NCBI Official Symbol	NOX4
NCBI Official Synonym Symbols	KOX; KOX-1; RENOX
NCBI Protein Information	NADPH oxidase 4
UniProt Protein Name	NADPH oxidase 4
UniProt Gene Name	NOX4
UniProt Synonym Gene Names	RENOX; KOX-1

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-rtu-elisa-kits/nicotinamide-adenine-dinucleotide-phosphate-oxidase-4-nox4/217426>