

Human Niemann Pick Disease Type C2 (NPC2) RTU ELISA Kit | NPC2 rtu elisa kit

Human Niemann Pick Disease Type C2 (NPC2) ELISA Kit

CATALOG #
AAA217429

REACTIVITY
Human


PRODUCT SYNONYMS

Niemann Pick Disease Type C2 (NPC2); Human Niemann Pick Disease Type C2 (NPC2) ELISA Kit; HE1; NP-C2; EDDM1; Epididymal Protein 1; Human epididymis-specific protein 1; Epididymal secretory protein E1; NPC2 rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	151
Assay Type	Sandwich
Samples	Serum, Plasma, Tissue Homogenates or Other Biological Fluids.
Detection Range	0.781-50ng/mL
Sensitivity	0.31ng/mL
RTU ELISA Kit's Advantage	The Ready-To-Use (RTU) Kits require fewer steps in setup and take less time to run. Detection Solutions are pre-diluted and total incubation time is shorter than many other ELISA Kits.
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for NPC2 rtu elisa kit

Test Principle: The microtiter plate provided in this kit has been pre-coated with an antibody specific to NPC2. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to NPC2. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain NPC2, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of NPC2 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	5453678
NCBI GeneID	10577
NCBI Accession #	NP_006423.1
NCBI GenBank Nucleotide #	NM_006432.4
UniProt Accession #	P61916; Q15668; Q29413; B4DQV7
Molecular Weight	13,416 Da
NCBI Official Full Name	NPC intracellular cholesterol transporter 2 isoform 2
NCBI Official Synonym Full Names	NPC intracellular cholesterol transporter 2
NCBI Official Symbol	NPC2
NCBI Official Synonym Symbols	HE1; EDDM1
NCBI Protein Information	NPC intracellular cholesterol transporter 2
UniProt Protein Name	NPC intracellular cholesterol transporter 2
UniProt Gene Name	NPC2
UniProt Synonym Gene Names	He1

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-rtu-elisa-kits/niemann-pick-disease-type-c2-npc2/217429>