

Canine Fbg ELISA Kit | Fbg elisa kit

Canine Fbg (Fibrinogen) ELISA Kit

CATALOG #
AAA21745

REACTIVITY
Canine


PRODUCT SYNONYMS

Fbg; Canine Fbg (Fibrinogen) ELISA Kit; Fbg elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Specificity	This kit recognizes Canine FG in samples. No significant cross-reactivity or interference between Canine FG and analogues was observed.
Sequence Length	128
Assay Type	Quantitative Sandwich
Samples	Serum, plasma and other biological fluids
Detection Range	62.5-4000pg/mL
Sensitivity	37.5pg/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, mid range and high level Canine FG were tested 20 times on one plate, respectively.
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, mid range and high level Canine FG were tested on 3 different plates, 20 replicates in each plate.
Preparation and Storage	Store at 4 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Fbg elisa kit	Principle of the Assay: This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Canine FG. Standards or samples are added to the micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Canine FG and Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Canine FG, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm ± 2 nm. The OD value is proportional to the concentration of Canine FG. You can calculate the concentration of Canine FG in the samples by comparing the OD of the samples to the standard curve. Intended Uses: This ELISA kit applies to the in vitro quantitative determination of Canine FG concentrations in serum, plasma and other biological fluids.
---	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	49247126
Molecular Weight	15,251 Da
NCBI Official Full Name	FREP4 FBG variant 3, partial
UniProt Protein Name	FREP4 FBG variant 3
UniProt Entry Name	Q5YD86_BIOGL

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/canine-elisa-kits/fbg/21745>