

Mouse Polyamine Oxidase (PAOX) RTU ELISA Kit | PAOX rtu elisa kit

Mouse Polyamine Oxidase (PAOX) ELISA Kit

CATALOG #
AAA217627

REACTIVITY
Mouse


PRODUCT SYNONYMS

Polyamine Oxidase (PAOX); Mouse Polyamine Oxidase (PAOX) ELISA Kit; PAO; Peroxisomal N(1)-Acetyl-Spermine/Spermidine Oxidase(Exo-N4-Amino); PAOX rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Sequence Length	504
Assay Type	Sandwich
Samples	Serum, Plasma, Tissue Homogenates or Other Biological Fluids.
Detection Range	78.1-5000pg/mL
Sensitivity	32pg/mL
RTU ELISA Kit's Advantage	The Ready-To-Use (RTU) Kits require fewer steps in setup and take less time to run. Detection Solutions are pre-diluted and total incubation time is shorter than many other ELISA Kits.
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PAOX rtu elisa kit

Test Principle: The microtiter plate provided in this kit has been pre-coated with an antibody specific to PAOX. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to PAOX. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain PAOX, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of PAOX in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	28173566
NCBI GeneID	212503
NCBI Accession #	NP_722478.2
NCBI GenBank Nucleotide #	NM_153783.4
UniProt Accession #	Q8C0L6; Q8K254
Molecular Weight	25,262 Da
NCBI Official Full Name	peroxisomal N(1)-acetyl-spermine/spermidine oxidase isoform 1
NCBI Official Synonym Full Names	polyamine oxidase (exo-N4-amino)
NCBI Official Symbol	Paox
NCBI Official Synonym Symbols	Pao; AI118225; 2410012F02Rik
NCBI Protein Information	peroxisomal N(1)-acetyl-spermine/spermidine oxidase
UniProt Protein Name	Peroxisomal N(1)-acetyl-spermine/spermidine oxidase
UniProt Gene Name	Paox
UniProt Synonym Gene Names	Pao

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/mouse-rtu-elisa-kits/polyamine-oxidase-paox/217627>