

Human Phosphatidylinositol Binding Clathrin Assembly Protein (PICALM) RTU ELISA Kit | PICALM rtu elisa kit

Human Phosphatidylinositol Binding Clathrin Assembly Protein (PICALM) ELISA Kit

CATALOG #
AAA217786

REACTIVITY
Human


PRODUCT SYNONYMS

Phosphatidylinositol Binding Clathrin Assembly Protein (PICALM); Human Phosphatidylinositol Binding Clathrin Assembly Protein (PICALM) ELISA Kit; LAP; CALM; CLTH; Clathrin assembly lymphoid myeloid leukemia protein; PICALM rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	610
Assay Type	Sandwich
Samples	Tissue Homogenates, Cell Lysates or Other Biological Fluids.
Detection Range	0.312-20ng/mL
Sensitivity	0.126ng/mL
RTU ELISA Kit's Advantage	The Ready-To-Use (RTU) Kits require fewer steps in setup and take less time to run. Detection Solutions are pre-diluted and total incubation time is shorter than many other ELISA Kits.
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PICALM rtu elisa kit

Test Principle: The microtiter plate provided in this kit has been pre-coated with an antibody specific to PICALM. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to PICALM. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain PICALM, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of PICALM in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	56788368
NCBI GeneID	8301
NCBI Accession #	NP_001008660.1
NCBI GenBank Nucleotide #	NM_001008660.2
UniProt Accession #	Q13492; O60700; Q4LE54; Q6GMQ6; Q86XZ9; B4DTM3; E9PN05; F8VPG7
Molecular Weight	70,037 Da
NCBI Official Full Name	phosphatidylinositol-binding clathrin assembly protein isoform 2
NCBI Official Synonym Full Names	phosphatidylinositol binding clathrin assembly protein
NCBI Official Symbol	PICALM
NCBI Official Synonym Symbols	LAP; CALM; CLTH
NCBI Protein Information	phosphatidylinositol-binding clathrin assembly protein
UniProt Protein Name	Phosphatidylinositol-binding clathrin assembly protein
UniProt Gene Name	PICALM
UniProt Synonym Gene Names	CALM

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-rtu-elisa-kits/phosphatidylinositol-binding-clathrin-assembly-protein-picalm/217786>