

Human IL-35 ELISA Kit | IL-35 elisa kit

Human IL-35 (Interleukin 35) ELISA Kit

CATALOG #
AAA21790

REACTIVITY
Human


PRODUCT SYNONYMS

IL-35; Human IL-35 (Interleukin 35) ELISA Kit; IL-35 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This kit recognizes Human IL-35 in samples. No significant cross-reactivity or interference between Human IL-35 and analogues was observed.
Sequence Length	229
Samples	Serum, plasma and other biological fluids
Assay Type	Quantitative Sandwich
Detection Range	15.63-1000pg/mL
Sensitivity	9.38pg/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, mid range and high level Human IL-35 were tested 20 times on one plate, respectively.
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, mid range and high level Human IL-35 were tested on 3 different plates, 20 replicates in each plate.
Preparation and Storage	Store at 4 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for IL-35 elisa kit

Intended Uses: This ELISA kit applies to the in vitro quantitative determination of Human IL-35 concentrations in serum, plasma and other biological fluids.

Principle of the Assay: This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Human IL-35. Standards or samples are added to the micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Human IL-35 and Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Human IL-35, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm +/- 2 nm. The OD value is proportional to the concentration of Human IL-35. You can calculate the concentration of Human IL-35 in the samples by comparing the OD of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	14577917
NCBI GeneID	10148
NCBI Accession #	NP_005746.2
NCBI GenBank Nucleotide #	NM_005755.2
UniProt Accession #	Q14213; O75269; A0N0N2
Molecular Weight	25,396 Da
NCBI Official Full Name	interleukin-27 subunit beta
NCBI Official Synonym Full Names	Epstein-Barr virus induced 3
NCBI Official Symbol	EBI3
NCBI Official Synonym Symbols	IL27B; IL-27B
NCBI Protein Information	interleukin-27 subunit beta; EBV-induced gene 3 protein; Epstein-Barr virus induced gene 3; IL-27 subunit beta; IL27 subunit; IL35 subunit; cytokine receptor; epstein-Barr virus-induced gene 3 protein
UniProt Protein Name	Interleukin-27 subunit beta
UniProt Gene Name	EBI3
UniProt Synonym Gene Names	IL27B; IL-27 subunit beta; IL-27B; EBV-induced gene 3 protein
UniProt Entry Name	IL27B_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/il-35/21790>