

Human Peroxiredoxin 5 (PRDX5) RTU ELISA Kit | PRDX5 rtu elisa kit

Human Peroxiredoxin 5 (PRDX5) ELISA Kit

CATALOG #
AAA217973

REACTIVITY
Human


PRODUCT SYNONYMS

Peroxiredoxin 5 (PRDX5); Human Peroxiredoxin 5 (PRDX5) ELISA Kit; PLP; PRDX6; ACR1; PMP20; PRXV; SBB10; Thioredoxin Peroxidase PMP20; Peroxisomal Antioxidant Enzyme; TPx Type VI; Liver Tissue 2D-Page Spot 71B; Alu Co-Repressor 1; PRDX5 rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	214
Assay Type	Sandwich
Samples	Tissue Homogenates, Cerebrospinal Fluid or Other Biological Fluids.
Detection Range	0.781-50ng/mL
Sensitivity	0.34ng/mL
RTU ELISA Kit's Advantage	The Ready-To-Use (RTU) Kits require fewer steps in setup and take less time to run. Detection Solutions are pre-diluted and total incubation time is shorter than many other ELISA Kits.
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PRDX5 rtu elisa kit

Test Principle: The microtiter plate provided in this kit has been pre-coated with an antibody specific to PRDX5. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to PRDX5. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain PRDX5, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of PRDX5 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	6912238
NCBI GeneID	25824
NCBI Accession #	NP_036226.1
NCBI GenBank Nucleotide #	NM_012094.4
UniProt Accession #	P30044; Q14CK0; Q6IAF2; Q9UBU5; Q9UJU4; Q9UKX4; A6NC19; A6NG06; B7ZLJ4; B7ZVW3
Molecular Weight	12,897 Da
NCBI Official Full Name	peroxiredoxin-5, mitochondrial isoform L
NCBI Official Synonym Full Names	peroxiredoxin 5
NCBI Official Symbol	PRDX5
NCBI Official Synonym Symbols	PLP; ACR1; B166; PRXV; PMP20; PRDX6; prx-V; SBBI10; AOEB166; HEL-S-55
NCBI Protein Information	peroxiredoxin-5, mitochondrial
UniProt Protein Name	Peroxiredoxin-5, mitochondrial
UniProt Gene Name	PRDX5
UniProt Synonym Gene Names	ACR1; AOEB166; Prx-V

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-rtu-elisa-kits/peroxiredoxin-5-prdx5/217973>