

Human Proteasome Activator Subunit 3 (PSME3) RTU ELISA Kit | PSME3 rtu elisa kit

Human Proteasome Activator Subunit 3 (PSME3) ELISA Kit

CATALOG #
AAA218032

REACTIVITY
Human


PRODUCT SYNONYMS

Proteasome Activator Subunit 3 (PSME3); Human Proteasome Activator Subunit 3 (PSME3) ELISA Kit; Ki; PA28G; REG-GAMMA; 11S regulator complex subunit gamma; Activator of multicatalytic protease subunit 3; Ki nuclear autoantigen; Proteasome activator 28 subunit gamma; PSME3 rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	254
Assay Type	Sandwich
Samples	Tissue Homogenates, Cell Lysates or Other Biological Fluids.
Detection Range	0.156-10ng/mL
Sensitivity	0.054ng/mL
RTU ELISA Kit's Advantage	The Ready-To-Use (RTU) Kits require fewer steps in setup and take less time to run. Detection Solutions are pre-diluted and total incubation time is shorter than many other ELISA Kits.
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PSME3 rtu elisa kit

Test Principle: The microtiter plate provided in this kit has been pre-coated with an antibody specific to PSME3. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to PSME3. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain PSME3, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of PSME3 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	30410794
NCBI GeneID	10197
NCBI Accession #	NP_005780.2
NCBI GenBank Nucleotide #	NM_005789.3
UniProt Accession #	P61289; O35563; P97373; Q12920; Q13172; Q9BQD9; A8K9A3
Molecular Weight	30,890 Da
NCBI Official Full Name	proteasome activator complex subunit 3 isoform 1
NCBI Official Synonym Full Names	proteasome activator subunit 3
NCBI Official Symbol	PSME3
NCBI Official Synonym Symbols	Ki; PA28G; HEL-S-283; PA28gamma; REG-GAMMA; PA28-gamma
NCBI Protein Information	proteasome activator complex subunit 3
UniProt Protein Name	Proteasome activator complex subunit 3
UniProt Gene Name	PSME3
UniProt Synonym Gene Names	REG-gamma; PA28g; PA28gamma

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-rtu-elisa-kits/proteasome-activator-subunit-3-psme3/218032>