

Human Retinoic Acid Receptor Responder 1 (RARRES1) RTU ELISA Kit | RARRES1 rtu elisa kit

Human Retinoic Acid Receptor Responder 1 (RARRES1) ELISA Kit

CATALOG #
AAA218131

REACTIVITY
Human


PRODUCT SYNONYMS

Retinoic Acid Receptor Responder 1 (RARRES1); Human Retinoic Acid Receptor Responder 1 (RARRES1) ELISA Kit; TIG1; PEIG1; Tazarotene-Induced Gene-1; Phorbol ester-induced gene 1 protein; RAR-responsive protein TIG1; RARRES1 rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	228
Assay Type	Sandwich
Samples	Serum, Plasma, Tissue Homogenates or Other Biological Fluids.
Detection Range	0.156-10ng/mL
Sensitivity	0.066ng/mL
RTU ELISA Kit's Advantage	The Ready-To-Use (RTU) Kits require fewer steps in setup and take less time to run. Detection Solutions are pre-diluted and total incubation time is shorter than many other ELISA Kits.
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for RARRES1 rtu elisa kit

Test Principle: The microtiter plate provided in this kit has been pre-coated with an antibody specific to RARRES1. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to RARRES1. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain RARRES1, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of RARRES1 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #

[46255041](#)

NCBI GeneID	5918
NCBI Accession #	NP_002879.2
NCBI GenBank Nucleotide #	NM_002888.3
UniProt Accession #	P49788; Q8N1D7
Molecular Weight	25,812 Da
NCBI Official Full Name	retinoic acid receptor responder protein 1 isoform 2
NCBI Official Synonym Full Names	retinoic acid receptor responder 1
NCBI Official Symbol	RARRES1
NCBI Official Synonym Symbols	LXNL; TIG1; PERG-1
NCBI Protein Information	retinoic acid receptor responder protein 1
UniProt Protein Name	Retinoic acid receptor responder protein 1
UniProt Gene Name	RARRES1
UniProt Synonym Gene Names	PEIG1; TIG1; PERG-1

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-rtu-elisa-kits/retinoic-acid-receptor-responder-1-rarres1/218131>