

Dog Surfactant Associated Protein A (SPA) RTU ELISA Kit | SPA rtu elisa kit

Canine Surfactant Associated Protein A (SPA) ELISA Kit

CATALOG #
AAA218527

REACTIVITY
Dog


PRODUCT SYNONYMS

Surfactant Associated Protein A (SPA); Canine Surfactant Associated Protein A (SPA) ELISA Kit; PRL; SFTPA; PSAP; PSPA; SFTP; SFTPA1; COLEC4; PSP-A; SFTP1; SFTPA1B; SP-A1; Collectin-4; Alveolar proteinosis protein; 35 kDa pulmonary surfactant-associated protein; SPA rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Dog
Sequence Length	262
Assay Type	Sandwich
Samples	Serum, Plasma, Lung Lavage Fluid, Tissue Homogenates or Other Biological Fluids.
Detection Range	1.87-120ng/mL
Sensitivity	0.62ng/mL
RTU ELISA Kit's Advantage	The Ready-To-Use (RTU) Kits require fewer steps in setup and take less time to run. Detection Solutions are pre-diluted and total incubation time is shorter than many other ELISA Kits.
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for SPA rtu elisa kit

Test Principle: The microtiter plate provided in this kit has been pre-coated with an antibody specific to SPA. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to SPA. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain SPA, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of SPA in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	346421445
NCBI GeneID	407213
NCBI Accession #	NP_001071306.2
NCBI GenBank Nucleotide #	NM_001077838.2
UniProt Accession #	Q6RXL1; Q08DV8
Molecular Weight	26,381 Da
NCBI Official Full Name	pulmonary surfactant-associated protein A
NCBI Official Symbol	SFTPA1
NCBI Official Synonym Symbols	PSAP; SP-A; PSP-A; SFTPA1B
NCBI Protein Information	pulmonary surfactant-associated protein A
UniProt Protein Name	Pulmonary surfactant-associated protein A
UniProt Gene Name	SFTPA1
UniProt Synonym Gene Names	SFTPA; PSAP; PSP-A; SP-A

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/dog-rtu-elisa-kits/surfactant-associated-protein-a-spa/218527>