

Human Serine Peptidase Inhibitor Kazal Type 5 (SPINK5) RTU ELISA Kit | SPINK5 rtu elisa kit

Human Serine Peptidase Inhibitor Kazal Type 5 (SPINK5) ELISA Kit

CATALOG # AAA218544	REACTIVITY Human
-------------------------------	----------------------------


PRODUCT SYNONYMS

Serine Peptidase Inhibitor Kazal Type 5 (SPINK5); Human Serine Peptidase Inhibitor Kazal Type 5 (SPINK5) ELISA Kit; LEKTI; LETKI; NETS; NS; VAKTI; Lymphoepithelial Kazal-Type-Related Inhibitor; Lympho-epithelial Kazal-type-related inhibitor; SPINK5 rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	1094
Assay Type	Sandwich
Samples	Tissue Homogenates, Cell Lysates or Other Biological Fluids.
Detection Range	0.312-20ng/mL
Sensitivity	0.116ng/mL
RTU ELISA Kit's Advantage	The Ready-To-Use (RTU) Kits require fewer steps in setup and take less time to run. Detection Solutions are pre-diluted and total incubation time is shorter than many other ELISA Kits.
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for SPINK5 rtu elisa kit	Test Principle: The microtiter plate provided in this kit has been pre-coated with an antibody specific to SPINK5. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to SPINK5. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain SPINK5, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of SPINK5 in the samples is then determined by comparing the O.D. of the samples to the standard curve.
--	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	189163501
-----------	-----------

NCBI GeneID	11005
NCBI Accession #	NP_001121170.1
NCBI GenBank Nucleotide #	NM_001127698.1
UniProt Accession #	Q9NQ38; O75770; Q3LX95; Q3LX96; Q3LX97; Q96PP2; Q96PP3; A8MYE8; B7WPB7; D6REN5
Molecular Weight	124,076 Da
NCBI Official Full Name	serine protease inhibitor Kazal-type 5 isoform a preproprotein
NCBI Official Synonym Full Names	serine peptidase inhibitor, Kazal type 5
NCBI Official Symbol	SPINK5
NCBI Official Synonym Symbols	NS; NETS; LEKTI; LETKI; VAKTI
NCBI Protein Information	serine protease inhibitor Kazal-type 5
UniProt Protein Name	Serine protease inhibitor Kazal-type 5
UniProt Gene Name	SPINK5
UniProt Synonym Gene Names	LEKTI

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-rtu-elisa-kits/serine-peptidase-inhibitor-kazal-type-5-spink5/218544>