

## Human Vitamin D Receptor (VDR) RTU ELISA Kit | VDR rtu elisa kit

Human Vitamin D Receptor (VDR) ELISA Kit

CATALOG #  
**AAA219154**

REACTIVITY  
**Human**


### PRODUCT SYNONYMS

Vitamin D Receptor (VDR); Human Vitamin D Receptor (VDR) ELISA Kit; NR111; Nuclear Receptor Subfamily 1 Group I Member 1; Calcitriol Receptor; 1; 25-dihydroxyvitamin D3 receptor; VDR rtu elisa kit

### PRODUCT SPECIFICATIONS

|                           |                                                                                                                                                                                                                                                                                                                                                     |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Reactivity                | Human                                                                                                                                                                                                                                                                                                                                               |
| Sequence Length           | 427                                                                                                                                                                                                                                                                                                                                                 |
| Assay Type                | Sandwich                                                                                                                                                                                                                                                                                                                                            |
| Samples                   | Tissue Homogenates, Cell Lysates or Other Biological Fluids.                                                                                                                                                                                                                                                                                        |
| Detection Range           | 0.625-40ng/mL                                                                                                                                                                                                                                                                                                                                       |
| Sensitivity               | 0.224ng/mL                                                                                                                                                                                                                                                                                                                                          |
| RTU ELISA Kit's Advantage | The Ready-To-Use (RTU) Kits require fewer steps in setup and take less time to run. Detection Solutions are pre-diluted and total incubation time is shorter than many other ELISA Kits.                                                                                                                                                            |
| Preparation and Storage   | For unopened kits: All the reagents should be kept at 4 degree C upon receipt.<br>For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal. |

### IMAGES



Standard Curve (Sample)

### ADDITIONAL INFORMATION

#### Related Product Information for VDR rtu elisa kit

**Test Principle:** The microtiter plate provided in this kit has been pre-coated with an antibody specific to VDR. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to VDR. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain VDR, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm  $\pm$  10nm. The concentration of VDR in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

|                                  |                                  |
|----------------------------------|----------------------------------|
| NCBI GI #                        | 4507883                          |
| NCBI GeneID                      | 7421                             |
| NCBI Accession #                 | NP_000367.1                      |
| NCBI GenBank Nucleotide #        | NM_000376.2                      |
| UniProt Accession #              | P11473; Q5PSV3; B2R5Q1; G3V1V9   |
| Molecular Weight                 | 53,883 Da                        |
| NCBI Official Full Name          | vitamin D3 receptor isoform VDRA |
| NCBI Official Synonym Full Names | vitamin D receptor               |
| NCBI Official Symbol             | VDR                              |
| NCBI Official Synonym Symbols    | NR1I1; PPP1R163                  |
| NCBI Protein Information         | vitamin D3 receptor              |
| UniProt Protein Name             | Vitamin D3 receptor              |
| UniProt Gene Name                | VDR                              |
| UniProt Synonym Gene Names       | NR1I1; VDR                       |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)  
Online datasheet: <https://stage.aaabiotech.com/human-rtu-elisa-kits/vitamin-d-receptor-vdr/219154>