

Human Xanthine Dehydrogenase (XDH) RTU ELISA Kit | XDH rtu elisa kit

Human Xanthine Dehydrogenase (XDH) ELISA Kit

CATALOG #
AAA219268

REACTIVITY
Human


PRODUCT SYNONYMS

Xanthine Dehydrogenase (XDH); Human Xanthine Dehydrogenase (XDH) ELISA Kit; XO; XOR; XOD; Xanthine Oxidase; NAD⁺-Xanthine Dehydrogenase; Xanthine-NAD⁺ Oxidoreductase; Xanthine/NAD⁺ Oxidoreductase; Xanthine Oxidoreductase; XDH rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	1333
Assay Type	Sandwich
Samples	Serum, Plasma, Breast Milk or Other Biological Fluids.
Detection Range	0.156-10ng/mL
Sensitivity	0.055ng/mL
RTU ELISA Kit's Advantage	The Ready-To-Use (RTU) Kits require fewer steps in setup and take less time to run. Detection Solutions are pre-diluted and total incubation time is shorter than many other ELISA Kits.
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for XDH rtu elisa kit

Test Principle: The microtiter plate provided in this kit has been pre-coated with an antibody specific to XDH. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to XDH. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain XDH, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of XDH in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	91823271
NCBI GeneID	7498
NCBI Accession #	NP_000370.2
NCBI GenBank Nucleotide #	NM_000379.3
UniProt Accession #	P47989; Q16681; Q16712; Q4PJ16
Molecular Weight	146,424 Da
NCBI Official Full Name	xanthine dehydrogenase/oxidase
NCBI Official Synonym Full Names	xanthine dehydrogenase
NCBI Official Symbol	XDH
NCBI Official Synonym Symbols	XO; XOR; XAN1
NCBI Protein Information	xanthine dehydrogenase/oxidase
UniProt Protein Name	Xanthine dehydrogenase/oxidase
UniProt Gene Name	XDH
UniProt Synonym Gene Names	XDHA; XO; XOR

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-rtu-elisa-kits/xanthine-dehydrogenase-xdh/219268>