

Rabbit anti-Human PERK Polyclonal Antibody | anti-PERK antibody

Rabbit anti-PERK Antibody

CATALOG #	HOST	REACTIVITY	APPLICATIONS
AAA223195	Rabbit	Human	Immunoprecipitation, Immunofluorescence, Immunohistochemistry, Western Blot


PRODUCT SYNONYMS

PERK; Polyclonal Antibody; Rabbit anti-PERK Antibody; PEK; Eukaryotic translation initiation factor 2-alpha kinase 3; PRKR-like endoplasmic reticulum kinase; Pancreatic eIF2-alpha kinase; HsPEK; anti-PERK antibody

PRODUCT SPECIFICATIONS


Host	Rabbit
Reactivity	Human
Clonality	Polyclonal
Purity/Purification	The antibody was purified by immunogen affinity chromatography.
Form/Format	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Concentration	1mg/mL (varies by lot)
Applicable Applications for anti-PERK antibody	IP (Immunoprecipitation), IF (Immunofluorescence), IHC (Immunohistochemistry), WB (Western Blot)
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human PERK. The exact sequence is proprietary.
Preparation and Storage	Store at -20 degree C. Avoid repeated freeze/thaw cycles.

IMAGES




IP (Immunoprecipitation)

Immunoprecipitation of PERK from 0.5mg HeLa whole cell extract lysate, using 5ug of anti-PERK Antibody and 50ul of protein G magnetic beads (+). No antibody was added to the control (-). The antibody was incubated under agitation with Protein G beads for 10min, HeLa whole cell extract lysate diluted in RIPA buffer was added to each sample and incubated for a further 10min under agitation. Proteins were eluted by addition of 40ul SDS loading buffer and incubated for 10min at 70°C; 10ul of each sample was separated on a SDS PAGE gel, transferred to a nitrocellulose membrane, blocked with 5% BSA and probed with anti-PERK Antibody.




IF (Immunofluorescence)

Immunofluorescent analysis of PERK staining in HeLa cells. Formalin-fixed cells were permeabilized with 0.1% Triton X-100 in TBS for 5-10 minutes and blocked with 3% BSA-PBS for 30 minutes at room temperature. Cells were probed with the primary antibody in 3% BSA-PBS and incubated overnight at 4 °C in a humidified chamber. Cells were washed with PBST and incubated with a DyLight 594-conjugated secondary antibody (red) in PBS at room temperature in the dark. DAPI was used to stain the cell nuclei (blue).



IHC (Immunohistochemistry)

Immunohistochemical analysis of PERK staining in human brain formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.



WB (Western Blot)

Western blot analysis of PERK expression in HeLa (A) whole cell lysates.

ADDITIONAL INFORMATION

Product Categories/Family for anti-PERK antibody

Primary antibodies

NCBI AND UNIPROT INFORMATION

NCBI GI #	134304838
NCBI GeneID	9451
NCBI Accession #	NP_004827.4
NCBI GenBank Nucleotide #	NM_004836.5
UniProt Accession #	Q9NZJ5; O95846; Q53QY0; Q53SB1; A0AVH1; A0AVH2; B2RCU9
Molecular Weight	125,216 Da
NCBI Official Full Name	eukaryotic translation initiation factor 2-alpha kinase 3
NCBI Official Synonym Full Names	eukaryotic translation initiation factor 2-alpha kinase 3
NCBI Official Symbol	EIF2AK3
NCBI Official Synonym Symbols	PEK; WRS; PERK
NCBI Protein Information	eukaryotic translation initiation factor 2-alpha kinase 3; hsPEK; pancreatic EIF2-alpha kinase; PRKR-like endoplasmic reticulum kinase; eukaryotic translation initiation factor 2 alpha kinase 3
UniProt Protein Name	Eukaryotic translation initiation factor 2-alpha kinase 3
UniProt Gene Name	EIF2AK3
UniProt Synonym Gene Names	PEK; PERK; HsPEK
UniProt Entry Name	E2AK3_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/polyclonal/perk-antibody-immunoprecipitation-immunofluorescence-immunohistochemistry-western-blot-elisa/223195>