

Rabbit p56Dok2 Polyclonal Antibody | anti-p56Dok2 antibody

p56Dok2 antibody (Tyr299)

CATALOG # AAA224478	HOST Rabbit	APPLICATIONS Immunohistochemistry, Western Blot
-------------------------------	-----------------------	---

PRODUCT SYNONYMS

p56Dok2; Polyclonal Antibody; p56Dok2 antibody (Tyr299); Polyclonal p56Dok2; Anti-p56Dok2; Docking protein 2; docking protein 2 56kDa; DOK2; Downstream of tyrosine kinase 2; p56(dok-2); p56DOK; p56dok-2; anti-p56Dok2 antibody

PRODUCT SPECIFICATIONS

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Specificity	p56Dok2 antibody (Tyr299) detects endogenous levels of p56Dok2 only when phosphorylated at tyrosine 299
Purity/Purification	p56Dok2 antibody was purified by affinity chromatography
Form/Format	Supplied in PBS buffer (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol
Concentration	1 mg/ml (varies by lot)
Applicable Applications for anti-p56Dok2 antibody	IHC (Immunohistochemistry), WB (Western Blot)
Biological Significance	p56Dok2 is a member of the p62dok family of proteins that interact with receptor tyrosine kinases and mediate particular biological responses.
Immunogen	p56Dok2 antibody was raised in Rabbit using a synthesized phosphopeptide derived from human p56Dok2 around the phosphorylation site of tyrosine 299 (G-E-YP-A-V) as the immunogen.
Preparation and Storage	Store at -20 degree C in small aliquots. Avoid repeat freeze/thaw cycles

IMAGES



WB (Western Blot)

Western blot analysis of extracts from K562 cells, using p56Dok-2 antibody and p56Dok2 antibody (Tyr299) (Lanes 3 and 4)



IHC (Immunohistochemistry)

Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using p56Dok2 antibody (Tyr299)

ADDITIONAL INFORMATION

Related Product Information for anti-p56Dok2 antibody	Rabbit Polyclonal p56Dok2 antibody (Tyr299)
Product Categories/Family for anti-p56Dok2 antibody	Cell Biology; Purified Polyclonal Antibodies

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

