

Canine Leptin ELISA Kit | Lep elisa kit

Canine Leptin (Lep) ELISA Kit

CATALOG #
AAA22568

REACTIVITY
Canine


PRODUCT SYNONYMS

Leptin; Canine Leptin (Lep) ELISA Kit; Lep elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Specificity	No cross-reaction with other factors.
Sequence Length	167
Assay Type	Sandwich
Detection Range	10 ng/mL-0.156 ng/mL
Sensitivity	Up to 0.05 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Lep elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Canine Lep monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	4557715
NCBI GeneID	3952

NCBI Accession #	NP_000221.1
NCBI GenBank Nucleotide #	NM_000230.2
Molecular Weight	18,641 Da
NCBI Official Full Name	leptin
NCBI Official Synonym Full Names	leptin
NCBI Official Symbol	LEP
NCBI Official Synonym Symbols	OB; OBS; LEPD
NCBI Protein Information	leptin
UniProt Protein Name	Leptin
UniProt Gene Name	LEP
UniProt Entry Name	LEP_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/lep-canine-elisa-kits/leptin/22568>