

Canine Testosterone ELISA Kit | HSD17B3 elisa kit

Canine Testosterone ELISA Kit

CATALOG #
AAA22589

REACTIVITY
Canine


PRODUCT SYNONYMS

Testosterone; Canine Testosterone ELISA Kit ; HSD17B3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Sequence Length	310
Assay Type	Sandwich
Detection Range	10 ng/mL-0.156 ng/mL
Sensitivity	Up to 0.05 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for HSD17B3 elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Canine Testo monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	4557649
NCBI GeneID	3293

NCBI Accession #	NP_000188.1
NCBI GenBank Nucleotide #	NM_000197.1
Molecular Weight	29,109 Da
NCBI Official Full Name	testosterone 17-beta-dehydrogenase 3
NCBI Official Synonym Full Names	hydroxysteroid (17-beta) dehydrogenase 3
NCBI Official Symbol	HSD17B3
NCBI Official Synonym Symbols	EDH17B3; SDR12C2
NCBI Protein Information	testosterone 17-beta-dehydrogenase 3
UniProt Protein Name	Testosterone 17-beta-dehydrogenase 3
UniProt Gene Name	HSD17B3
UniProt Synonym Gene Names	EDH17B3; SDR12C2; 17-beta-HSD 3
UniProt Entry Name	DHB3_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/canine-elisa-kits/testosterone/22589>