

Canine Epidermal Growth Factor ELISA Kit | EGF elisa kit

Canine Epidermal Growth Factor (EGF) ELISA Kit

CATALOG #
AAA22616

REACTIVITY
Canine


PRODUCT SYNONYMS

Epidermal Growth Factor; Canine Epidermal Growth Factor (EGF) ELISA Kit; EGF elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Sequence Length	53
Assay Type	Sandwich
Detection Range	1000 pg/mL-15.6 pg/mL
Sensitivity	Up to 5 pg/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for EGF elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Canine EGF monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	374922206
NCBI GeneID	1950

Molecular Weight	129,746 Da
NCBI Official Full Name	epidermal growth factor, partial
NCBI Official Synonym Full Names	epidermal growth factor
NCBI Official Symbol	EGF
NCBI Official Synonym Symbols	URG; HOMG4
NCBI Protein Information	pro-epidermal growth factor
UniProt Protein Name	Pro-epidermal growth factor
UniProt Gene Name	EGF
UniProt Synonym Gene Names	EGF
UniProt Entry Name	EGF_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/egf-canine-elisa-kits/epidermal-growth-factor/22616>