

Mouse Pro Brain Natriuretic Peptide ELISA Kit | ProBNP elisa kit

Mouse Pro Brain Natriuretic Peptide (ProBNP) ELISA Kit

CATALOG #
AAA22668

REACTIVITY
Mouse


PRODUCT SYNONYMS

Pro Brain Natriuretic Peptide; Mouse Pro Brain Natriuretic Peptide (ProBNP) ELISA Kit; ProBNP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	No cross-reaction with other factors.
Sequence Length	760
Samples	Mouse serum, plasma or cell culture supernatant
Assay Type	Quantitative Sandwich
Detection Range	2000 pg/ml-31.2 pg/ml
Sensitivity	Up to 12 pg/ml
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for ProBNP elisa kit

Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is Mouse ProBNP monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	292495099
NCBI GeneID	2191

Molecular Weight	26,954 Da
NCBI Official Full Name	Prolyl endopeptidase FAP
NCBI Official Synonym Full Names	fibroblast activation protein alpha
NCBI Official Symbol	FAP
NCBI Official Synonym Symbols	FAPA; SIMP; DPPIV
NCBI Protein Information	prolyl endopeptidase FAP
UniProt Protein Name	Prolyl endopeptidase FAP
UniProt Gene Name	FAP
UniProt Synonym Gene Names	APCE
UniProt Entry Name	SEPR_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/probnp-mouse-elisa-kits/pro-brain-natriuretic-peptide/22668>