

Bovine Collagen Type III (COL3) ELISA Kit | COL3 elisa kit

Bovine Collagen Type III (COL3) ELISA Kit

CATALOG #
AAA22713

REACTIVITY
Bovine


PRODUCT SYNONYMS

Collagen Type III (COL3); Bovine Collagen Type III (COL3) ELISA Kit; COL3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Bovine
Sequence Length	1466
Assay Type	Sandwich
Detection Range	20 ng/mL-0.312 ng/mL
Sensitivity	Up to 0.06 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for COL3 elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The precoated antibody is an anti-Bovine COL3 monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidinperoxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	119631316
NCBI GeneID	1281
Molecular Weight	111,907 Da

NCBI Official Full Name	collagen, type III, alpha 1 (Ehlers-Danlos syndrome type IV, autosomal dominant), isoform CRA_b
NCBI Official Synonym Full Names	collagen type III alpha 1 chain
NCBI Official Symbol	COL3A1
NCBI Official Synonym Symbols	EDS4A
NCBI Protein Information	collagen alpha-1(III) chain
UniProt Protein Name	Collagen alpha-1(III) chain
UniProt Gene Name	COL3A1
UniProt Entry Name	CO3A1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/bovine-ovine-elisa-kits/collagentypeiii-col3/22713>