

Hemagglutinin glycoprotein Recombinant Protein | H recombinant protein

Recombinant Measles virus Hemagglutinin glycoprotein (H)

CATALOG # AAA230786	HOST Cell Free Expression
-------------------------------	-------------------------------------

PRODUCT SYNONYMS

Hemagglutinin glycoprotein; Recombinant Measles virus Hemagglutinin glycoprotein (H); H recombinant protein

PRODUCT SPECIFICATIONS

Host	Cell Free Expression
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Liquid containing glycerol
Sequence Positions	1-617. Full length.
Sequence	MSPQRDRINA FYKDNPHPKG SRVINREHL MIDRPYVLLA VLFVMFLSLI GLLAIAGIRL HRAAIYTAEI HKSLSTNLDV TNSIEHQVKD VLTPLFKIIG DEVGLRTPQR FTDLVKFISD KIKFLNPDRE YDFRDLTWCI NPPERIKLDY DQYCADVAE ELMNALVNST LLETRTTNQF LAVSKGNCSSG PTTIRGQFSN MSLSLDLYL GRGYNVSSIV TMTSQGMYGG TYLVEKPNLS SKRSELSQLS MYRVFEVGV I RNPGLGAPVF HMTNYLEQPV SNDLSNCMVA LGELKLAALC HGEDSITIPY QGSGKGVSFQ LVKLGWVKSP TDMQSWVPLS TDDPVIDRLY LSSHRGVIAD NQAKWAVPTT RTDDKLRLMET CFQQACKGKI QALCENPEWA PLKDNRIPSY GVLSVDLSLT VELKIKIASG FGPLITHGSG MDLYKSNHNN VYWLTIIPPMK NLALGVINTL EWIPRFKVSP YLFNVPIKEA GEDCHAPTYL PAEVDGDVKL SSNLVILPGQ DLQYVLATYD TSRVEHAVVY YVYSPSRFS YFYFPLPIK GVPIELQVEC FTWDQKLWCR HFCVLADSES GGHITHSGME GMGVSCTVTR EDGTNRR
Species	Measles virus (strain Edmonston) (MeV) (Subacute sclerose panencephalitis virus)
Preparation and Storage	Store at -20 degree C, for extended storage, conserve at -20 degree C or -80 degree C. Repeated freezing and thawing is not recommended. Store working aliquots at 4 degree C for up to one week.

ADDITIONAL INFORMATION

Product Categories/Family for H recombinant protein	Transmembrane Protein
---	---------------------------------------

NCBI AND UNIPROT INFORMATION

NCBI GI #	122985
UniProt Accession #	P08362
Molecular Weight	69,250 Da
NCBI Official Full Name	Hemagglutinin glycoprotein
UniProt Protein Name	Hemagglutinin glycoprotein
UniProt Gene Name	H
UniProt Entry Name	HEMA_MEASE

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

