

Schlafen family member 12 (SLFN12) Recombinant Protein | SLFN12 recombinant protein

Recombinant Human Schlafen family member 12 (SLFN12)

CATALOG # AAA233487	HOST Cell Free Expression
-------------------------------	-------------------------------------


PRODUCT SYNONYMS

Schlafen family member 12 (SLFN12); Recombinant Human Schlafen family member 12 (SLFN12); SLFN12 recombinant protein

PRODUCT SPECIFICATIONS

Host	Cell Free Expression
Purity/Purification	Greater or equal to 85% purity as determined by SDS-PAGE.
Form/Format	Liquid containing glycerol
Sequence Positions	1-578aa; Full Length Protein
Sequence	MNISVDLETN YAEVLVDVGR VTLGENSRKK MKDCKLRKKQ NESVSRAMCA LLNSGGGVK AEIENEDYSY TKDGIGLDLE NSFSNILLFV PEYLDFMQNG NYFLIFVKSX SLNTSGLRIT TLSSNLYKRD ITSAKVMNAT AALEFLKDMK KTRGRLYLRP ELLAKRPCVD IQEENNMKAL AGVFFDRTEL DRKEKLTFTS THVEIKNFS TEKLLQRIKE ILPQYVSAFA NTDGGYLFIG LNEDKEIIGF KAEMSDLDDL EREIEKSIRK MPVHHFCMEK KKINYSCKFL GYVDKGSCLG YVCALRVERF CCAVFAKEPD SWHVKDNRVM QLTRKEWQF MVEAEPKFSS SYEEVISQIN TSLPAPHSWP LLEWQRQRHH CPGLSGRITY TPENLCRKLQ LQHEGLKQLI CEEMDSVRKG SLIFSRWSV DLGLQENHKV LCDALLISQD SPPVLYTFHM VQDEEFGYS TQTALTLKQK LAKIGGYTKK VCVMTKIFYL SPEGMTSCQY DLRSQVIYPE SYFFTRRKYL LKALFKALKR LKSLRDQFSF AENLYQIIGI DCFQKNDKKM FKSCRRLT
Preparation and Storage	Store at -20 degree C, for extended storage, conserve at -20 degree C or -80 degree C. Repeated freezing and thawing is not recommended. Store working aliquots at 4 degree C for up to one week.

IMAGES



SDS-PAGE

ADDITIONAL INFORMATION

Related Product Information for SLFN12 recombinant protein	Interacting selectively and non-covalently with ATP, adenosine 5'-triphosphate, a universally important coenzyme and enzyme regulator.
Product Categories/Family for SLFN12 recombinant protein	Transmembrane Protein

NCBI AND UNIPROT INFORMATION

Molecular Weight	83.0 kDa
------------------	----------

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

