

Human Vasorin/SLIT-Like 2 ELISA Kit | VASN elisa kit

Human Vasorin/SLIT-Like 2 ELISA Kit

CATALOG #
AAA23618

REACTIVITY
Human


PRODUCT SYNONYMS

Vasorin/SLIT-Like 2; Human Vasorin/SLIT-Like 2 ELISA Kit; VASN elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	673
Sensitivity	1.0 pg/mL
Intra-assay Precision	Intra-assay CV is less than 15%.
Inter-assay Precision	Inter-assay CV is less than 15%.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for VASN elisa kit

Intended Uses: This VASORIN;SLIT-LIKE 2 ELISA kit is intended Laboratory for Research use only and is not for use in diagnostic or therapeutic procedures. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of VASORIN;SLIT-LIKE 2 in the sample, this VASORIN;SLIT-LIKE 2 ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus VASORIN;SLIT-LIKE 2 concentration. The concentration of VASORIN;SLIT-LIKE 2 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	88702793
NCBI GenBank Nucleotide #	NP_612449.2

UniProt Accession #	Q6EMK4
NCBI Official Full Name	vasorin
UniProt Protein Name	Vasorin
UniProt Gene Name	VASN
UniProt Synonym Gene Names	SLITL2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/vasorin-slit-like-2/23618>