

Human Heparin Sodium ELISA Kit

Human Heparin Sodium ELISA Kit

CATALOG #
AAA23662

REACTIVITY
Human


PRODUCT SYNONYMS

Heparin Sodium; Human Heparin Sodium ELISA Kit; Heparin Sodium elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Detection Range	7.5 ng/mL - 240 ng/mL
Sensitivity	1.0 ng/mL
Intra-assay Precision	Intra-assay CV are less than 15%.
Inter-assay Precision	Inter-assay CV are less than 15%.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Heparin Sodium elisa kit

Intended Uses: This HEPARIN SODIUM ELISA kit is intended Laboratory for Research use only and is not for use in diagnostic or therapeutic procedures. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of HEPARIN SODIUM in the sample, this HEPARIN SODIUM ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus HEPARIN SODIUM concentration. The concentration of HEPARIN SODIUM in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/heparin-sodium/23662>