

General Penicillin ELISA Kit

General Penicillin ELISA Kit

CATALOG #
AAA23664

REACTIVITY
General


PRODUCT SYNONYMS

Penicillin; General Penicillin ELISA Kit; Penicillin elisa kit

PRODUCT SPECIFICATIONS

Reactivity	General
Sensitivity	1 ng/ml
Intra-assay Precision	Intra-assay CV are less than 15%.
Inter-assay Precision	Inter-assay CV are less than 15%.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Penicillin elisa kit

Intended Uses: This PENICILLIN ELISA kit is intended Laboratory for Research use only and is not for use in diagnostic or therapeutic procedures. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of PENICILLIN in the sample, this PENICILLIN ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus PENICILLIN concentration. The concentration of PENICILLIN in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.