

Human Heparin-PF4 Complex; HIT Antibody ELISA Kit

Human Heparin-PF4 Complex; HIT Antibody ELISA Kit

CATALOG #
AAA23691

REACTIVITY
Human


PRODUCT SYNONYMS

Heparin-PF4 Complex; HIT Antibody; Human Heparin-PF4 Complex; HIT Antibody ELISA Kit; HIT Antibody elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sensitivity	0.1 ng/mL
Intra-assay Precision	Intra-assay CV is less than 15%.
Inter-assay Precision	Inter-assay CV is less than 15%.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Heparin-PF4 Complex; HIT Antibody elisa kit

Intended Uses: This HEPARIN-PF4 COMPLEX; HIT ANTIBODY ELISA kit is intended Laboratory for Research use only and is not for use in diagnostic or therapeutic procedures. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of HEPARIN-PF4 COMPLEX; HIT ANTIBODY in the sample, this HEPARIN-PF4 COMPLEX; HIT ANTIBODY ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus HEPARIN-PF4 COMPLEX; HIT ANTIBODY concentration. The concentration of HEPARIN-PF4 COMPLEX; HIT ANTIBODY in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.