

Human Salivary Oxytocin ELISA Kit

Human Salivary Oxytocin ELISA Kit

CATALOG #
AAA23700

REACTIVITY
Human


PRODUCT SYNONYMS

Salivary Oxytocin; Human Salivary Oxytocin ELISA Kit; Salivary Oxytocin elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Detection Range	25 pg/mL - 800 pg/mL
Sensitivity	1 pg/mL
Intra-assay Precision	Intra-assay CV are less than 15%.
Inter-assay Precision	Inter-assay CV are less than 15%.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Salivary Oxytocin elisa kit

Intended Uses: This SOT ELISA kit is intended Laboratory for Research use only and is not for use in diagnostic or therapeutic procedures. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. To measure the concentration of SOT in the sample, this SOT ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a standard curve of Optical Density versus SOT concentration. The concentration of SOT in the samples is then determined by comparing the O.D. of the samples to the standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/salivary-oxytocin/23700>