

Human Protein S (PROS) RTU ELISA Kit | PROS rtu elisa kit

Human Protein S (PROS) ELISA Kit

CATALOG #
AAA23979

REACTIVITY
Human


PRODUCT SYNONYMS

Protein S (PROS); Human Protein S (PROS) ELISA Kit; PSA; PS; PS26; PS21; PS22; PS23; PS24; PS25; VKDPS; VKDP-S; Protein Sa; Vitamin K Dependent Protein S; PROS rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of PROS. No significant cross-reactivity or interference between PROS and analogues was observed.
Sequence Length	676
Assay Type	Quantitative Sandwich
Samples	Human plasma
Detection Range	1.56-100ng/mL
Sensitivity	< 0.73ng/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level PROS were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level PROS were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PROS rtu elisa kit	Principle of the Assay: The microtiter plate provided in this kit has been pre-coated with an antibody specific to PROS. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to PROS. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain PROS, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of PROS in the samples is then determined by comparing the O.D. of the samples to the standard curve. Intended Uses: The kit is a sandwich enzyme immunoassay for the in vitro quantitative measurement of PROS in human plasma.
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	192447438
NCBI GeneID	5627
NCBI Accession #	NP_000304.2
NCBI GenBank Nucleotide #	NM_000313.3
UniProt Accession #	P07225; Q15518; Q7Z715; Q9UCZ8; A8KAC9; D3DN28
Molecular Weight	75,123 Da
NCBI Official Full Name	vitamin K-dependent protein S isoform 2 preproprotein
NCBI Official Synonym Full Names	protein S
NCBI Official Symbol	PROS1
NCBI Official Synonym Symbols	PSA; PROS; PS21; PS22; PS23; PS24; PS25; THPH5; THPH6
NCBI Protein Information	vitamin K-dependent protein S
UniProt Protein Name	Vitamin K-dependent protein S
UniProt Gene Name	PROS1
UniProt Synonym Gene Names	PROS

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.