

Rabbit degP Polyclonal Antibody | anti-degP antibody

Rabbit anti-Escherichia coli (strain K12) degP Polyclonal antibody

CATALOG # AAA244160	HOST Rabbit	REACTIVITY E Coli (strain K12)	APPLICATIONS Western Blot, ELISA
-------------------------------	-----------------------	--	--


PRODUCT SYNONYMS

degP; Polyclonal Antibody; Rabbit anti-Escherichia coli (strain K12) degP Polyclonal antibody; degP Antibody; Periplasmic serine endoprotease DegP (EC 3.4.21.107)(Heat shock protein DegP)(Protease Do); htrA ptd; anti-degP antibody

PRODUCT SPECIFICATIONS

Host	Rabbit
Reactivity	E Coli (strain K12)
Clonality	Polyclonal
Isotype	IgG
Purity/Purification	Protein A/G
Form/Format	Liquid Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH7.4
Applicable Applications for anti-degP antibody	WB (Western Blot), ELISA
Immunogen	Recombinant Escherichia coli Periplasmic serine endoprotease DegP (degP)(27-474AA)
Immunogen Species	Escherichia coli (strain K12)
Preparation and Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

IMAGES



WB (Western Blot)

Western Blot Positive WB detected in Recombinant protein All lanes: degP antibody at 1:2000Secondary Goat polyclonal to rabbit IgG at 1/50000 dilutionPredicted band size: 62.8 kDaObserved band size: 65 kDa

ADDITIONAL INFORMATION

Product Categories/Family for anti-degP antibody	Polyclonal Antibodies
--	-----------------------

NCBI AND UNIPROT INFORMATION

NCBI GI #	16128154
NCBI GeneID	947139
NCBI Accession #	NP_414703.1
NCBI GenBank Nucleotide #	NC_000913.3
UniProt Accession #	P0C0V0; P09376; P15724
NCBI Official Full Name	serine endoprotease (protease Do), membrane-associated
NCBI Official Symbol	degP
NCBI Official Synonym Symbols	ECK0160; htrA; JW0157; ptd
NCBI Protein Information	serine endoprotease (protease Do), membrane-associated
UniProt Protein Name	Periplasmic serine endoprotease DegP
UniProt Gene Name	degP
UniProt Synonym Gene Names	htrA; ptd
UniProt Entry Name	DEGP_ECOLI

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.