

Human anti SSA/Ro antibody ELISA Kit

Human anti SSA/Ro antibody ELISA Kit

CATALOG #
AAA244437

REACTIVITY
Human


PRODUCT SYNONYMS

anti SSA/Ro antibody; Human anti SSA/Ro antibody ELISA Kit; anti SSA/Ro antibody elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Detection Range	25-500ng/mL
Sensitivity	1.0ng/mL
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for anti SSA/Ro antibody elisa kit

Intended Uses: This SSA/Ro Ab ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Human SSA/Ro Ab. This ELISA kit for research use only!

Principle of the Assay: SSA/Ro Ab ELISA kit applies the competitive enzyme immunoassay technique utilizing a SSA/Ro antigen and an SSA/Ro Ab-HRP conjugate. The assay sample and buffer are incubated together with SSA/Ro Ab-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the SSA/Ro Ab concentration since SSA/Ro Ab from samples and SSA/Ro Ab-HRP conjugate compete for the SSA/Ro antigen binding site. Since the number of sites is limited, as more sites are occupied by SSA/Ro Ab from the sample, fewer sites are left to bind SSA/Ro Ab-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The SSA/Ro Ab concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/anti-ssa-ro-antibody/244437>