

Canine anti-cyclic citrullinated peptide antibody (anti-CCP antibody) ELISA Kit

Canine anti-cyclic citrullinated peptide antibody (anti-CCP antibody) ELISA Kit

CATALOG #
AAA244476

REACTIVITY
Canine


PRODUCT SYNONYMS

anti-cyclic citrullinated peptide antibody (anti-CCP antibody); Canine anti-cyclic citrullinated peptide antibody (anti-CCP antibody) ELISA Kit; anti-cyclic citrullinated peptide antibody (anti-CCP antibody) elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Samples	Serum, Plasma, Cell Culture Supernatants, Body Fluid And Tissue Homogenate
Assay Type	Quantitative Competitive
Sensitivity	1.0 pg/mL.
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for anti-cyclic citrullinated peptide antibody (anti-CCP antibody) elisa kit

Intended Uses: This CCL-26 ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Bovine CCL-26. This ELISA kit for research use only!

Principle of the Assay: CCL-26 ELISA kit applies the competitive enzyme immunoassay technique utilizing a polyclonal anti-CCL-26 antibody and an CCL-26-HRP conjugate. The assay sample and buffer are incubated together with CCL-26-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the CCL-26 concentration since CCL-26 from samples and CCL-26-HRP conjugate compete for the anti-CCL-26 antibody binding site. Since the number of sites is limited, as more sites are occupied by CCL-26 from the sample, fewer sites are left to bind CCL-26-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The CCL-26 concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.