

Human Glutaredoxin 1 (GLRX) ELISA Kit | GLRX elisa kit

Human Glutaredoxin 1 (GLRX) ELISA Kit

CATALOG #
AAA245297

REACTIVITY
Human


PRODUCT SYNONYMS

Glutaredoxin 1 (GLRX); Human Glutaredoxin 1 (GLRX) ELISA Kit; GLRX elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Detection Range	250-5000pg/mL
Sensitivity	1.0pg/mL
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for GLRX elisa kit

Intended Uses: This GLRX ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Human GLRX.

Principle of the Assay: GLRX ELISA kit applies the competitive enzyme immunoassay technique utilizing a Polyclonal anti-GLRX antibody and an GLRX-HRP conjugate. The assay sample and buffer are incubated together with GLRX-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the GLRX concentration since GLRX from samples and GLRX-HRP conjugate compete for the anti-GLRX antibody binding site. Since the number of sites is limited, as more sites are occupied by GLRX from the sample, fewer sites are left to bind GLRX-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The GLRX concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/glutaredoxin-1-glrx/245297>