

Canine Irisin ELISA Kit | Irisin elisa kit

Canine Irisin ELISA Kit

CATALOG #
AAA245426

REACTIVITY
Canine


PRODUCT SYNONYMS

Irisin; Canine Irisin ELISA Kit; Irisin elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Samples	Serum, Plasma, Cell Culture Supernatants, Body Fluid And Tissue Homogenate
Assay Type	Quantitative Competitive
Sensitivity	0.1 ng/mL.
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Irisin elisa kit

Intended Uses: This IRF3 ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Canine IRF3. This ELISA kit for research use only!

Principle of the Assay: IRF3 ELISA kit applies the competitive enzyme immunoassay technique utilizing an anti-IRF3 antibody and an IRF3-HRP conjugate. The assay sample and buffer are incubated together with IRF3-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the IRF3 concentration since IRF3 from samples and IRF3-HRP conjugate compete for the anti-IRF3 antibody binding site. Since the number of sites is limited, as more sites are occupied by IRF3 from the sample, fewer sites are left to bind IRF3-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The IRF3 concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/canine-elisa-kits/irisin/245426>