

Bovine Adherens junction-associated protein 1 (AJAP1) ELISA Kit | AJAP1 elisa kit

Bovine Adherens junction-associated protein 1 (AJAP1) ELISA Kit

CATALOG # AAA245965	REACTIVITY Bovine
-------------------------------	-----------------------------


PRODUCT SYNONYMS

Adherens junction-associated protein 1 (AJAP1); Bovine Adherens junction-associated protein 1 (AJAP1) ELISA Kit; AJAP1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Bovine
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for AJAP1 elisa kit	Intended Uses: This AIRA ELISA kit is intended Laboratory for Research use only. The Stop Solution changes the color from blue to yellow and the intensity of the color is measured at 450 nm using a spectrophotometer. In order to measure the concentration of AIRA in the sample, this AIRA ELISA Kit includes a set of calibration standards. The calibration standards are assayed at the same time as the samples and allow the operator to produce a cutoff value. The existence or not of AIRA in the samples is then determined by comparing the O.D. of the samples to the CUT OFF.
Product Categories/Family for AJAP1 elisa kit	Cancer

NCBI AND UNIPROT INFORMATION

NCBI GI #	109715839
NCBI GeneID	55966
NCBI Accession #	NP_001035943.1
NCBI GenBank Nucleotide #	NM_001042478.1
Molecular Weight	44,536 Da
NCBI Official Full Name	adherens junction-associated protein 1
NCBI Official Synonym Full Names	adherens junctions associated protein 1
NCBI Official Symbol	AJAP1
NCBI Official Synonym Symbols	MOT8; SHREW1; SHREW-1

NCBI Protein Information	adherens junction-associated protein 1; membrane protein shrew-1; transmembrane protein SHREW1; adherens junction associated protein 1
UniProt Protein Name	Adherens junction-associated protein 1
UniProt Gene Name	AJAP1
UniProt Synonym Gene Names	MOT8; SHREW1
UniProt Entry Name	AJAP1_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/bovine-ovine-elisa-kits/adherens-junction-associated-protein-1-ajap1/245965>