

Canine Endomorphin-2 ELISA Kit | EM-2 elisa kit

Canine Endomorphin-2 ELISA Kit

CATALOG #
AAA246084

REACTIVITY
Canine


PRODUCT SYNONYMS

Endomorphin-2; Canine Endomorphin-2 ELISA Kit; EM-2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Samples	Serum, plasma, cell culture supernatants, body fluid and tissue homogenate
Assay Type	Quantitative Competitive
Sensitivity	1.0 pg/mL.
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



ADDITIONAL INFORMATION

Related Product Information for EM-2 elisa kit

Intended Uses: This EM2 ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Canine EM2. This ELISA kit for research use only!

Principle of the Assay: EM2 ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-EM2 antibody and an EM2-HRP conjugate. The assay sample and buffer are incubated together with EM2-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the EM2 concentration since EM2 from samples and EM2-HRP conjugate compete for the anti-EM2 antibody binding site. Since the number of sites is limited, as more sites are occupied by EM2 from the sample, fewer sites are left to bind EM2-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The EM2 concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

