

Mouse Adenosine Monophosphate ELISA Kit | AMP elisa kit

Mouse Adenosine Monophosphate ELISA Kit

CATALOG #
AAA246265

REACTIVITY
Mouse


PRODUCT SYNONYMS

Adenosine Monophosphate; Mouse Adenosine Monophosphate ELISA Kit; AMP elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	This assay has high sensitivity and excellent specificity for detection of AMP. No significant cross-reactivity or interference between AMP and analogues was observed. NOTE: Limited by current skills and knowledge, it is impossible for us to complete the cross-reactivity detection between AMP and all the analogues, therefore, cross reaction may still exist in some cases.
Assay Type	Quantitative Competitive
Samples	Serum, plasma, cell culture supernatants, body fluid and tissue homogenate
Sensitivity	1.0 ng/mL
Intended Uses	This AMP ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Mouse AMP. This ELISA kit for research use only, not for therapeutic applications!
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for AMP elisa kit	Principle of the Assay: AMP ELISA kit applies the competitive enzyme immunoassay technique utilizing a polyclonal anti-AMP antibody and an AMP-HRP conjugate. The assay sample and buffer are incubated together with AMP-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the AMP concentration since AMP from samples and AMP-HRP conjugate compete for the anti-AMP antibody binding site. Since the number of sites is limited, as more sites are occupied by AMP from the sample, fewer sites are left to bind AMP-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The AMP concentration in each sample is interpolated from this standard curve. Intended Uses: This AMP ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Mouse AMP. This ELISA kit for research use only, not for therapeutic or test applications!
Product Categories/Family for AMP elisa kit	Immunology

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/amp-mouse-elisa-kits/adenosine-monophosphate/246265>