

Canine CD117 ELISA Kit | CD117 elisa kit

Canine CD117 ELISA Kit

CATALOG #
AAA246312

REACTIVITY
Canine


PRODUCT SYNONYMS

CD117; Canine CD117 ELISA Kit; CD117 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Sensitivity	0.1 ng/mL.
Intended Uses	This CD117 ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Canine CD117. This ELISA kit for research use only, not for therapeutic applications!
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for CD117 elisa kit

Principle of the assay: CD117 ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-CD117 antibody and an CD117-HRP conjugate. The assay sample and buffer are incubated together with CD117-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the CD117 concentration since CD117 from samples and CD117-HRP conjugate compete for the anti-CD117 antibody binding site. Since the number of sites is limited, as more sites are occupied by CD117 from the sample, fewer sites are left to bind CD117-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The CD117 concentration in each sample is interpolated from this standard curve.

Product Categories/Family for CD117 elisa kit

Immunology

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

