

Human blca 4 ELISA Kit

Human blca 4 ELISA Kit

CATALOG #
AAA246512

REACTIVITY
Human


PRODUCT SYNONYMS

blca 4; Human blca 4 ELISA Kit; blca 4 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, cell culture supernatants, body fluid and tissue homogenate
Assay Type	Competitive (Quantitative)
Sensitivity	1.0 mg/L.
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for blca 4 elisa kit

Intended Uses: This UTM ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Human UTM. This ELISA kit for research use only!

Principle of the Assay: UTM ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-UTM antibody and an UTM-HRP conjugate. The assay sample and buffer are incubated together with UTM-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the UTM concentration since UTM from samples and UTM-HRP conjugate compete for the anti-UTM antibody binding site. Since the number of sites is limited, as more sites are occupied by UTM from the sample, fewer sites are left to bind UTM-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The UTM concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.