

Human Tumor Necrosis Factor Receptor 2 ELISA Kit | TNFR2 elisa kit

Human Tumor Necrosis Factor Receptor 2 ELISA Kit

CATALOG #
AAA246992

REACTIVITY
Human


PRODUCT SYNONYMS

Tumor Necrosis Factor Receptor 2; Human Tumor Necrosis Factor Receptor 2 ELISA Kit; TNFR2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Serum, plasma, Cell Culture Supernatants, body fluid and tissue homogenate
Assay Type	Competitive
Detection Range	250-5000pg/mL
Sensitivity	1.0pg/mL
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for TNFR2 elisa kit

Intended Uses: This TNFR2 ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Human TNFR2. This ELISA kit for research use only!

Principle of the Assay: TNFR2 ELISA kit applies the competitive enzyme immunoassay technique utilizing a polyclonal anti-TNFR2 antibody and an TNFR2-HRP conjugate. The assay sample and buffer are incubated together with TNFR2-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the TNFR2 concentration since TNFR2 from samples and TNFR2-HRP conjugate compete for the anti-TNFR2 antibody binding site. Since the number of sites is limited, as more sites are occupied by TNFR2 from the sample, fewer sites are left to bind TNFR2-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The TNFR2 concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://stage.aaabiotech.com/tnfr2-human-elisa-kits/tumor-necrosis-factor-receptor-2/246992>