

Canine Urinary Retinol Binding Protein ELISA Kit

Canine Urinary Retinol Binding Protein ELISA Kit

CATALOG #
AAA247472

REACTIVITY
Canine


PRODUCT SYNONYMS

Urinary Retinol Binding Protein; Canine Urinary Retinol Binding Protein ELISA Kit; Urinary Retinol Binding Protein elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Canine
Specificity	This assay has high sensitivity and excellent specificity for detection of RBP. No significant cross-reactivity or interference between RBP and analogues was observed. NOTE: Limited by current skills and knowledge, it is impossible for us to complete the cross-reactivity detection between RBP and all the analogues, therefore, cross reaction may still exist in some cases.
Samples	Serum, plasma, cell culture supernatants, body fluid and tissue homogenate
Assay Type	Quantitative Competitive
Sensitivity	1.0 ng/mL
Reactivity Note	No significant cross-reactivity or interference between this antigen and analogues was observed. NOTE: Limited by current skills and knowledge, it is impossible for us to complete the cross-reactivity detection, therefore, cross reaction may still exist in some cases.
Preparation and Storage	Store at 2-8 degree C. The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Urinary Retinol Binding Protein elisa kit	Intended Uses: This RBP ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Canine RBP. This ELISA kit for research use only! Principle of the Assay: RBP ELISA kit applies the competitive enzyme immunoassay technique utilizing an anti-RBP antibody and an RBP-HRP conjugate. The assay sample and buffer are incubated together with RBP-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the RBP concentration since RBP from samples and RBP-HRP conjugate compete for the anti-RBP antibody binding site. Since the number of sites is limited, as more sites are occupied by RBP from the sample, fewer sites are left to bind RBP-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The RBP concentration in each sample is interpolated from this standard curve.
--	--

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/canine-elisa-kits/urinary-retinol-binding-protein/247472>