

Human Anti-cochlear antibody ELISA Kit

Human Anti-cochlear antibody ELISA Kit

CATALOG #
AAA247593

REACTIVITY
Human


PRODUCT SYNONYMS

Anti-cochlear antibody; Human Anti-cochlear antibody ELISA Kit; Anti-cochlear antibody elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of ACAb. No significant cross-reactivity or interference between ACAb and analogues was observed. NOTE: Limited by current skills and knowledge, it is impossible for us to complete the cross-reactivity detection between ACAb and all the analogues, therefore, cross reaction may still exist in some cases.
Assay Type	Quantitative Competitive
Samples	Serum, plasma, cell culture supernatants, body fluid and tissue homogenate
Sensitivity	1.0 ng/mL
Intra-Assay	CV<10%
Inter-Assay	CV<15%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for Anti-cochlear antibody elisa kit	Principle of the Assay: ACAb ELISA kit applies the competitive enzyme immunoassay technique utilizing Cochlear antigen and an ACAb-HRP conjugate. The assay sample and buffer are incubated together with ACAb-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the ACAb concentration since ACAb from samples and ACAb-HRP conjugate compete for the Cochlear antigen binding site. Since the number of sites is limited, as more sites are occupied by ACAb from the sample, fewer sites are left to bind ACAb-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The ACAb concentration in each sample is interpolated from this standard curve. Intended Uses: This ACAb ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Human ACAb. This ELISA kit for research use only, not for therapeutic or test applications!
---	---

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/anti-cochlear-antibody/247593>