

Bovine Guanidinoacetate N-methyltransferase ELISA Kit | GAMT elisa kit

Bovine Guanidinoacetate N-methyltransferase ELISA Kit

CATALOG # AAA247854	REACTIVITY Bovine	APPLICATIONS ELISA
-------------------------------	-----------------------------	------------------------------


PRODUCT SYNONYMS

Guanidinoacetate N-methyltransferase; Bovine Guanidinoacetate N-methyltransferase ELISA Kit; GAMT elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Bovine
Specificity	This assay has high sensitivity and excellent specificity for detection of GAMT. No significant cross-reactivity or interference between GAMT and analogues was observed. NOTE: Limited by current skills and knowledge, it is impossible for us to complete the cross-reactivity detection between GAMT and all the analogues, therefore, cross reaction may still exist in some cases.
Assay Type	Quantitative Competitive
Samples	Serum, plasma, cell culture supernatants, body fluid and tissue homogenate
Sensitivity	1.0 ng/mL
Applicable Applications for GAMT elisa kit	ELISA
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for GAMT elisa kit	Principle of the Assay: GAMT ELISA kit applies the competitive enzyme immunoassay technique utilizing an anti-GAMT antibody and an GAMT-HRP conjugate. The assay sample and buffer are incubated together with GAMT-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the GAMT concentration since GAMT from samples and GAMT-HRP conjugate compete for the anti-GAMT antibody binding site. Since the number of sites is limited, as more sites are occupied by GAMT from the sample, fewer sites are left to bind GAMT-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The GAMT concentration in each sample is interpolated from this standard curve. Intended Uses: This GAMT ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Bovine GAMT. This ELISA kit for research use only, not for therapeutic or test applications!
--	---

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/gamt-bovine-ovine-elisa-kits/guanidinoacetate-n-methyltransferase/247854>