

Human Anti-bullous pemphigoid 230 antibody IgE (BP230-Ab IgE) ELISA Kit | BP230-Ab IgE elisa kit

Human Anti-bullous pemphigoid 230 antibody IgE (BP230-Ab IgE) ELISA Kit

CATALOG #
AAA247918

REACTIVITY
Human


PRODUCT SYNONYMS

Anti-bullous pemphigoid 230 antibody IgE (BP230-Ab IgE); Human Anti-bullous pemphigoid 230 antibody IgE (BP230-Ab IgE) ELISA Kit; BP230-Ab IgE elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of BP230-Ab IgE. No significant cross-reactivity or interference between BP230-Ab IgE and analogues was observed. NOTE: Limited by current skills and knowledge, it is impossible for us to complete the cross-reactivity detection between BP230-Ab IgE and all the analogues, therefore, cross reaction may still exist in some cases.
Assay Type	Quantitative Competitive
Samples	Serum, plasma, cell culture supernatants, body fluid and tissue homogenate
Sensitivity	0.1 mg/ml
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for BP230-Ab IgE elisa kit	Principle of the Assay: BP230-Ab IgE ELISA kit applies the competitive enzyme immunoassay technique utilizing an anti-BP230-Ab IgE antibody and an BP230-Ab IgE-HRP conjugate. The assay sample and buffer are incubated together with BP230-Ab IgE-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the BP230-Ab IgE concentration since BP230-Ab IgE from samples and BP230-Ab IgE-HRP conjugate compete for the anti-BP230-Ab IgE antibody binding site. Since the number of sites is limited, as more sites are occupied by BP230-Ab IgE from the sample, fewer sites are left to bind BP230-Ab IgE-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The BP230-Ab IgE concentration in each sample is interpolated from this standard curve. Intended Uses: This BP230-Ab IgE ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Human BP230-Ab IgE. This ELISA kit for research use only, not for therapeutic or test applications!
---	--

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/anti-bullous-pemphigoid-230-antibody-ige-bp230-ab-ige/247918>