

Human INT6 (Integrator complex subunit 6) ELISA Kit | INT6 elisa kit

Human INT6 (Integrator complex subunit 6) ELISA Kit

CATALOG #
AAA248608

REACTIVITY
Human


PRODUCT SYNONYMS

INT6 (Integrator complex subunit 6); Human INT6 (Integrator complex subunit 6) ELISA Kit; Integrator complex subunit 6; INT6 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	854
Assay Type	Sandwich ELISA, Double Antibody
Samples	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Detection Range	78.125-5000pg/ml
Sensitivity	46.875pg/ml
Preparation and Storage	Store entire kit at 2-8C for short-term. For longer-term, please store the microplate & standard at -20C, while the remaining reagents can be stored at 2-8C

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for INT6 elisa kit

Principle of the Assay: This kit was based on sandwich enzyme-linked immune-sorbent assay technology. Anti-INT6 antibody was pre-coated onto 96-well plates. And the biotin conjugated anti-INT6 antibody was used as detection antibodies. The standards, test samples and biotin conjugated detection antibody were added to the wells subsequently, and wash with wash buffer. HRP-Streptavidin was added and unbound conjugates were washed away with wash buffer. TMB substrates were used to visualize HRP enzymatic reaction. TMB was catalyzed by HRP to produce a blue color product that changed into yellow after adding acidic stop solution. The density of yellow is proportional to the INT6 amount of sample captured in plate. Read the O.D. absorbance at 450nm in a microplate reader, and then the concentration of INT6 can be calculated.

NCBI AND UNIPROT INFORMATION

NCBI GI #

187608071

NCBI GeneID	569456
NCBI Accession #	NP_001119900.1
NCBI GenBank Nucleotide #	NM_001126428.1
NCBI Official Full Name	integrator complex subunit 6
NCBI Official Synonym Full Names	integrator complex subunit 6
NCBI Official Symbol	ints6
NCBI Official Synonym Symbols	ddx26a
NCBI Protein Information	integrator complex subunit 6

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/int6-integrator-complex-subunit-6/248608>