

Human MAGEA3 (Melanoma-associated antigen 3) ELISA Kit | MAGEA3 elisa kit

Human MAGEA3 (Melanoma-associated antigen 3) ELISA Kit

CATALOG #
AAA248842

REACTIVITY
Human

PRODUCT SYNONYMS

MAGEA3 (Melanoma-associated antigen 3); Human MAGEA3 (Melanoma-associated antigen 3) ELISA Kit; MAGEA3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of MAGEA3. No significant cross-reactivity or interference between MAGEA3 and analogues was observed.
Sequence Length	314
Assay Type	Sandwich ELISA, Double Antibody
Samples	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Detection Range	0.313-20ng/ml
Sensitivity	0.188ng/ml
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level MAGEA3 were tested 20 times on one plate, respectively. Intra-Assay: CV<8%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level MAGEA3 were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<10%
Preparation and Storage	Store entire kit at 2-8C for short-term. For longer-term, please store the microplate & standard at -20C, while the remaining reagents can be stored at 2-8C

IMAGES

Typical Data & Standard Curve

Results of a typical standard run of a MAGEA3 ELISA Kit are shown below. This standard curve was generated at our lab for demonstration purpose only. Each user should obtain their own standard curve as per experiment. (N/A: not applicable)

X	ng/ml	0	0.312	0.625	1.25	2.5	5	10	20
Y	OD450	0.024	0.088	0.152	0.257	0.489	0.887	1.331	1.919

Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for MAGEA3 elisa kit

Principle of the Assay: This kit was based on sandwich enzyme-linked immune-sorbent assay technology. Anti-MAGEA3 antibody was pre-coated onto 96-well plates. And the biotin conjugated anti-MAGEA3 antibody was used as detection antibodies. The standards, test samples and biotin conjugated detection antibody were added to the wells subsequently, and wash with wash buffer. HRP-Streptavidin was added and unbound conjugates were washed away with wash buffer. TMB substrates were used to visualize HRP enzymatic reaction. TMB was catalyzed by HRP to produce a blue color product that changed into yellow after adding acidic stop solution. The density of yellow is proportional to the MAGEA3 amount of sample captured in plate. Read the O.D. absorbance at 450nm in a microplate reader, and then the concentration of MAGEA3 can be calculated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	4885467
NCBI GeneID	4102
NCBI Accession #	NP_005353.1
NCBI GenBank Nucleotide #	NM_005362.3
UniProt Accession #	P43357; Q6FHI6
Molecular Weight	34,747 Da
NCBI Official Full Name	melanoma-associated antigen 3
NCBI Official Synonym Full Names	MAGE family member A3
NCBI Official Symbol	MAGEA3
NCBI Official Synonym Symbols	HIP8; HYPD; CT1.3; MAGE3; MAGEA6
NCBI Protein Information	melanoma-associated antigen 3
UniProt Protein Name	Melanoma-associated antigen 3
UniProt Gene Name	MAGEA3
UniProt Synonym Gene Names	MAGE3; CT1.3

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://stage.aaabiotech.com/human-elisa-kits/magea3-melanoma-associated-antigen-3/248842>